320 XT

320 XB

320 IP65

Service Manual – Serviceanleitung

Identification – Identifikation

Service manual 320 XT/XB/IP65
350-8100
Index d1
February 2012
Subject to technical modifications.

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www.precisa.com

Copyright

No reproduction in any form of this manual, in whole or in part, may be made without written authorization from Precisa Gravimetrics AG.

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© Precisa Gravimetrics AG, 8953 Dietikon, Switzerland, 2012

1.Mechanics – Mechanik	3
1.1 Balance groups and spare parts – Waagengruppen und Ersatzteile	4
1.2 Housing's bottom – Gehäuseunterteil	6
1.3 Housing's top – Gehäuseoberteil	8
1.4 Housing's top, IP65 – Gehäuseoberteil, IP65	11
1.5 Weighing cell – Wägezelle	12
1.6 Weighing pan D80 Group A2 – Waagschale D80 Gruppe A2	15
1.7 Weighing pan 135, Group M – Waagschale 135, Gruppe M	16
1.8 Weighing pan 170, Group C – Waagschale 170, Gruppe C	17
1.9 Weighing pan 200, Group D – Waagschale 200, Gruppe D	18
1.10 Windshield, XT – Windschutz, XT 350-8517	19
1.11 Windshield 180, XB – Windschutz 180, XB 350-8518	21
1.12 Windshield 260, XB – Windschutz 260, XB 350-8519	23
1.13 Windshield easy – Windschutz einfach 320-8504	25
1.14 Twin-Reference Weight Device – Zwei-Referenzgewichts-Modul 320-7140	26
1.15 Single Reference Weight Device – Ein-Referenzgewichts-Modul 320-8021	28
2.Electronics – Elektronik	30
2.1 Main board, new vers. – Hauptprint, neue Ausf. 320-7215-020/021	31
2.2 Main board, old version – Hauptprint, alte Ausführung 320-7200-020/021	31
2.3 Connector board – Steckerprint 320-7232-010	32
2.4 Connector board – Steckerprint 320-7203-010	32
2.5 Display board VFD – Displayprint VFD 320-7204-010	33
2.6 Display board LCD – Displayprint LCD 320-7207-010	33
2.7 Universal switching power adapter – Universal-Schaltnetzteil 350-8384	34
2.8 Universal switching power adapter IP65 – Universal-Schaltnetzteil IP65 320-8025	34
3.Twin Reference Weight Version	35
3.1 Before you get started	36
3.2 Opening a balance	39
3.3 Removing the weighing cell from the housing's bottom	40
3.4 Dismantling the weighing cell	41
3.5 Installing new horizontal flexures	43
3.6 Cleaning	44
3.7 Assembling the weighing cell	45
3.8 Installing the weighing cell into the housing's bottom	50
4.Single Reference Weight Version	51
4.1 Before you get started	52
4.2 Opening a balance	55
4.3 Removing the reference weight device (SCS only)	56
4.4 Removing the weighing cell from the housing's bottom	56
4.5 Dismantling the weighing cell	57
4.6 Installing new horizontal flexures	59
4.7 Cleaning	60
4.8 Assembling the weighing cell	61
4.9 Installing the weighing cell into the housing's bottom	66
4.10 Installing the reference weight device into the housing's bottom (SCS only)	66
5.Adjustments	67
5.1 Before you get started	68
5.2 Setting the balance to the factory mode (the balance is open)	68
5.3 Adjusting the symmetry	69

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Certain spare parts (e.g. flexures) have to be determined more precisely because of the various capacities of balances. Therefore refer to the following tables and find out to which group a balance belongs.

Example: A 620M belongs to group M1. Thus, in group specific spare parts select the ones of group M1.

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



All other parts see see "Housing's top – Gehäuseoberteil" on page 8

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
1.5.2 Flexure repair sets – Lager-Reparatursets

A repair set consists of the required spare flexures for one balance, i.e. 1 coupling, 2 vertical flexures and 8 horizontal flexures.

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
1.6 Weighing pan D80 Group A2 – Waagschale D80 Gruppe A2

1.6.1 Weighing pan, old version – Waagschale, alte Version

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

1.6.2 Weighing pan, new version – Waagschale, neue Version

For XT balances from series 2439 and all 28xx, for XB balances from series 3139 –
Für XT-Waagen ab Serie 2439 und alle 28xx, für XB-Waagen ab Serie 3139


Pos.	Article number	Article	Artikel
2	320-8040	Weighing pan, for windsh. XT	Waagschale kpl., für Windschutz XT
	320-8041	Weighing pan, for windsh. XB	Waagschale kpl., für Windschutz XB

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










2.1 Main board, new vers. – Hauptprint, neue Ausf. 320-7215-020/021

2.2 Main board, old version – Hauptprint, alte Ausführung 320-7200-020/021


The old basic version of the main board 320-7200-020 or 320-7200-021 is no more available. In case of need, use the new, fully compatible version. In case of special software, please contact our PRECISA customer services dept.

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


3.1 Before you get started


3.1.1 Introduction

1. In case of a broken flexure, all flexures must be replaced. Furthermore it is safest to put all parts in their pristine place when re-assembling, particularly the screws!
2. Enough space and a sorted way of working are best conditions.

3.1.2 Conventions these instructions use to follow

1. (26/p.12): component items with page reference, in round brackets.
[5]: tooling items as listed in the kit, see "Tool kit 350-8537" on page 37, in square brackets.
Example: 3 nuts (26/p.12) means 3 nuts as item 26 on page 12.

3.1.3 Tool kit 350-8537

Pos.	Article number	Article	Quantity
1	W 41-1688	Mounting plate	1
2	W 41-1688-5	Support post long	2
3	W 41-1688-6	Support post short	2
4	350-6203	Screwdriver Torx M4	1
5	240-7124	Screwdriver small No.1	1
6	PW 29.43.300	Pin for flexure holders	2
7	W 41-1700	Assembly jig with slot	1
8	W 41-1701	Assembly jig with hole	1
9	W 41-1677	Assembly pin for assembly jig	4
10	W 41-1691	Centring spigot for flexures	2
11	W 41-1774	Centring spigot for coil, EP/ES	1
12	W 41-1675	Centring spigot for coil, XT/XB	1
13	W 41-1674	Coupling anchor	1
14	PN 1100-172	Screw cyl. Torx M4x12	2
15	W 41-1699	Alignment jig	1
16	W 41-1578-8	Connecting cable	1
17	W 41-1681	Centring bush	1

3.1.4 Further tools needed


Instruments:

1. Soldering iron
2. Tweezers
3. Lens
4. Voltmeter

Standard tools:

1. Box spanner Nos.4, 5.5 and 7
2. Flat-bladed screwdriver Nos.1, 2 and 4
3. Phillips screwdriver No.2
4. Hexagon socket wrench WS 2.5mm

3.1.5 Preparing the tools


1. Screw the 2 long [2] and 2 short [3] support posts on the assembly plate.
2. Preheat the soldering iron to $310 \pm 5^{\circ}\text{C}$ ($590 \pm 10^{\circ}\text{F}$).

3.2 Opening a balance

3.2.1 General procedure for opening an XT/XB balance


1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 cylindrical screws on the housing's top using a No. 4 screwdriver.
4. Slightly raise the housing's top.
5. Release the membrane keypad ribbon cable from the main board.

3.2.2 General procedure for opening an IP65 balance

1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 screws IP56 with a No. 1 screwdriver and remove the clamping ring.
4. Unscrew the 3 cylindrical screws on the housing's top and the 2 cylindrical screws on the housing's bottom using a No. 4 screwdriver.
5. Slightly raise the housing's top (do not damage the sealing).
6. Release the membrane keypad ribbon cable from the main board.

3.3 Removing the weighing cell from the housing's bottom

1. Unplug the floating sensor (15/p.12) from the main board using a small screwdriver.
2. **SCS only:** Remove the 4 screws (12/p.26) on the reference weight device from the counterholder (4/p.26) and remove the counterholder and the reference weight.
3. Turn over the housing's bottom and unscrew the 3 nuts (26/p.12) on the underside of the housing
4. Hold the weighing cell only on the chassis and lift it out of the housing's bottom.
5. The chassis can now be placed on the assembly plate for easier handling.
6. To protect the weigh cell, fix the 2 assembly jigs [7] + [8] with the 4 assembly pins [9] .

3.4 Dismantling the weighing cell

3.4.1 Removing the coupling

1. Unscrew the 2 screws (25/p.12), remove the coupling (3/p.12) incl. its bases (9/p.12) and spring washers (28/p.12) → Do not separate the screw components.
2. Screw the coupling anchor [13] tightly on with the 2 screws [14] in place of the coupling.

3.4.2 Removing the 2 flexure holders

1. Fix the upper flexure holder (2/p.12) into the chassis, using the 2 pins for flexure holders [6].
2. Unscrew the 4 screws (24/p.12) of the upper flexure holder and remove them together with its 4 spring washers (28/p.12) → 2 on the chassis, 2 on the support piece.
3. Remove the 2 pins with the upper flexure holder together with its 4 spring washers and flexures, then remove the pins from the holder.
4. Fix the lower flexure holder (2/p.12) into the chassis with the two pins for flexure holders [6].
5. Unscrew the 4 screws (24/p.12) from the lower flexure holder and remove them together with the 4 spring washers → 2 on the chassis, 2 on the support piece.
6. Remove the 2 pins with the lower flexure holder together with its 4 spring washers and flexures, then remove the pins from the holder.

3.4.3 Removing the coil

1. Unsolder the 2 connecting wires (1/p.12) carefully from the connecting board (17/p.12) using the tweezers.
2. Unscrew the 3 screws (23/p.12) from the magnet cover (8/p.12), remove the cover together with the transport safety device (4/p.12) and the spring washer (28/p.12).
3. Unscrew the 2 brass screws (21/p.12) from the coil (12/p.12) on the balance arm (19/p.12) using a No.2 screwdriver, remove the coil together with the washer (27/p.12).
4. Remove the coil carefully from the pot (34/p.12).

3.4.4 Removing the support piece

1. Unscrew the 3 screws (23/p.12) of the inner cone (12/p.12) and remove it.
2. Remove the left (10/p.12) and the right (9/p.12) weight holder (XT only).
3. Unscrew the upper screw (balance arm side) of the coupling anchor [13] and pull off the support piece together with the coupling anchor over the assembly pins [9].
4. Untighten the lower screw of the coupling anchor for roughly 2 turns. → The support piece and the coupling anchor are now loosely held together.

3.4.5 Removing the balance arm


1. Unscrew the brass screw (30/p.12) from the stopper (29/p.12), turn the stopper inwards ca. 90° and gently re-tighten the brass screw.
2. Unscrew the 4 screws (25/p.12) from the vertical flexures (6/p.12) and remove them together with its bases (11/p.12) and spring washers (28/p.12). **Note:** Since the lug of the floating sensor is firmly attached to the balance arm, the floating sensor must be unscrewed to dismantle the balance arm and gently tilted back:
3. Unscrew the 2 screws (20/p.12) from the floating sensor plate (5/p.12) and remove the plate.
4. Unscrew the screw (20/p.12) from the connection board (17/p.12).
5. Unscrew the 3 screws (22/p.12) of the floating sensor and gently tilt the floating sensor to the rear.
6. Remove the two inner assembly pins [9] → unscrew from the balance arm.
7. Remove the two outer assembly pins [9] → unscrew from the chassis.
8. Remove the two assembly jigs [7] + [8].
9. Carefully tilt the balance arm (19/p.12) over the pot and the floating sensor and take it out.
10. Further dismantling of the weighing cell is unnecessary.

3.5 Installing new horizontal flexures

3.5.1 Removing the old horizontal flexures

1. Remove all screws from the flexure holders and remove the old flexures together with their spring washers.


3.5.2 Installing new horizontal flexures


1. Place a flexure holder on the assembly plate [1] .
2. Select 4 new horizontal flexures (7/p.12) → Place them carefully on the assembly plate → their slot has to be always on the flexure holder.
3. Place the 4 spring washers (28/p.12) on the flexures and screw them tightly to the holder with 4 screws (22/p.12).
4. Take the complete flexure holder carefully off the assembly plate → Do not touch the flexures!
5. Prepare the other flexure holder the same way.

3.6 Cleaning


3.6.1 Cleaning the pot

1. Wrap a strips a double-sided adhesive tape around a flat stick.

2. Dab with the stick on the ground of the pot.

3. Prepare a thin and round stick the same way and spin it through the gap between the pot's wall and the magnet.
4. Dab over the top of the magnet pot with a piece of adhesive tape.

3.6.2 Cleaning the coil


1. Using a lint-free cloth, remove any contamination from the coil → Do not clean the coil with adhesive tape; risk of injury.

3.6.3 Checking the cleanliness again

1. Inspect the annular gap and the coil for residual contamination and clean if necessary.

3.7 Assembling the weighing cell

3.7.1 Assembling the balance arm

1. To facilitate the assembly, the chassis can be placed on the assembly plate → Pay attention to the floating sensor!
2. Place the assembly jigs [7] + [8] with 2 assembly pins [9] on the chassis, position them horizontally using the alignment jig [15] and tighten the two assembly pins [9] with a No. 4 box spanner.
3. Carefully introduce the balance arm (19/p.12) over the pot (34/p.12) and past it to the floating sensor.
4. Assemble the balance arm with 2 assembly pins [9] loosely behind the assembly jigs. The balance arm must remain moveable.
5. The right inner assembly pin of the assembly jig [8] can now be firmly tightened.
6. The left inner assembly pin of the assembly jig [7] can now be firmly tightened.


3.7.2 Assembling the coil, the magnet cover and the transport safety device

1. Place the coil onto the balance arm, centre it with the centring spigot [12] and fix it tightly with 2 brass screws (21/p.12) and 2 washers (27/p.12).
2. The centring spigot must be smoothly removable.
3. Pre-fix the magnet cover (8/p.12) with its left and right screw (23/p.12) .
4. Pre-assemble the transport safety device (4/p.12) with the third screw (23/p.12) and spring washer (28/p.12) → tighten gently.
5. The centring bush [17] shall be introduced between the balance arm-pilot and the transport safety device to assure constant annular gap. Tighten the screw (23/p.12) and remove the centring bush.

3.7.3 Assembling the floating sensor and the stopper


1. Press the floating sensor carefully into the cutout of the chassis in the direction of the pot and secure it to the chassis with the 3 screws (22/p.12).
2. Screw firmly the connection board (17/p.12) onto the chassis with the screw (20/p.12).
3. Fix the floating sensor plate (5/p.12) with the 2 screws (20/p.12) to the chassis.
4. Solder the two connection wires (1/p.12) parallelly and strain-freely to the connection board (17/p.12).
5. Turn back the stopper (29/p.12) with its brass screw (30/p.12) by 90° and tighten it.
6. The stopper has to be free of play and lie centrally in the symmetry screw (32/p.12) → if necessary, adjust the symmetry screw.

3.7.4 Assembling new vertical flexures


1. Select 2 new vertical flexures (6/p.12) see "Balance groups and spare parts – Waagengruppen und Ersatzteile" on page 4.
2. Bring the flexures into position → its slots have to be on the support piece side (below).
3. Centre the flexures in the hole on the chassis side (above) with a centring spigot [10].
4. Centre the flexures in the slot on the support piece side (below) with the second centring spigot [10].
5. Tighten the first (upper) centring spigot and then remove the second (lower) centring spigot. → Do not touch the flexures!
6. Fasten the flexures on the balance arm side with their bases (11/p.12), spring washers (28/p.12) and screws (25/p.12).
7. Remove the second (upper) centring spigot and fasten the vertical flexures on the chassis side with their bases (11/p.12), spring washers (28/p.12) and screws (25/p.12).

3.7.5 Assembling the support piece, the inner cone and the weight holders


1. Introduce the support piece with the loosely held coupling anchor over the two inner assembly pins [9].
2. It must be possible to push the support piece smoothly over the two inner assembly pins. In case the support piece can not be pushed smoothly running over the two assembly pins, unscrew the pin [9] and tighten it again.
3. Locate the support piece and the balance arm using the coupling anchor [13] and the two screws [14] → The support piece and balance arm are fixed now.
4. Assemble the inner cone (12/p.12) to the support piece with the 3 screws (23/p.12).
5. **SCS only:** Assemble the two weight holders (9/p.12) and (10/p.12) .

3.7.6 Assembling the flexure holders

1. Introduce the two pins for flexure holders [6] in the first pre-assembled holder (→ see "Installing new horizontal flexures" on page 43) the way that the shoulder of the pins is opposite the screw heads on the holder.
2. Bring the two pins with the holder into position in the chassis.
3. Lower the holder carefully over the two pins on the chassis, until the flexures just touch the chassis and the support piece.
4. Fix the holder with the 4 screws (24/p.12) and the 4 spring washers (28/p.12).
5. Push the two pins through, up to the installed horizontal flexures holder.
6. Introduce on the opposite side the second pre-assembled flexure holder (→ see "Installing new horizontal flexures" on page 43) with the screw heads inwards over the two pins and lower onto the chassis, until the flexures just touch the chassis and the support piece.
7. Secure the holder with the 4 screws (24/p.12) and the 4 spring washers (28/p.12).
8. It must be possible to pull out the two pins smoothly from the flexure holders.

3.7.7 Assembling a new coupling

1. Select a new coupling (3/p.12) see "Balance groups and spare parts – Waagengruppen und Ersatzteile" on page 4.
2. Remove the coupling anchor [13].
3. Bring the coupling into position → the slot is on the support piece side.
4. Centre the coupling first on the balance arm side, then on the support piece side, using a centring spigot [10].
5. Tighten the first (upper) centring spigot, then carefully remove the second (lower) centring spigot. Do not touch the coupling and fasten it on the support piece side with the base (11/p.12), the spring washer (28/p.12) and the screw (25/p.12) .
6. Remove the second (upper) centring spigot and fasten the coupling on the balance arm side with the base (11/p.12), the spring washer (28/p.12) and the screw (25/p.12).

3.7.8 Dismantling the assembly jigs

1. Remove the 2 inner assembly pins [9] carefully.
2. Remove the 2 outer assembly pins [9] carefully.
3. Take out the two assembly jigs [7] + [8] by sliding them sideways.
4. Remove the two assembly pins for flexure holders [6] carefully.

3.8 Installing the weighing cell into the housing's bottom

1. Carry out a further movement check with the lens on the balance arm → the balance arm must not touch the transport safety device (4/p.12).
2. Place the weighing cell carefully into the housing's bottom → hold the weighing cell only on the chassis.
3. Tilt the housing's bottom to the side and screw the weighing cell securely from below with the 3 nuts (26/p.12).
4. Connect the floating sensor plug (15/p.12) to the main board (11/p.6).
5. **SCS only:** Insert the internal reference weights and mount the counterholders (4/p.26) with the 4 screws (12/p.26).

4. Single Reference Weight Version

4.1 Before you get started


4.1.1 Introduction

1. In case of a broken flexure, all flexures must be replaced. Furthermore it is safest to put all parts in their pristine place when re-assembling, particularly the screws!
2. Enough space and a sorted way of working are best conditions.

4.1.2 Conventions these instructions use to follow

1. (26/p.12): component items with page reference, in round brackets.
2. [5]: tooling items as listed in the kit, see "Tool kit 350-8537" on page 53, in square brackets.
Example: 3 nuts (26/p.12) means 3 nuts as item 26 on page 12.

4.1.3 Tool kit 350-8537

Pos.	Article number	Description	Quantity
1	W 41-1688	Mounting plate	1
2	W 41-1688-5	Support post long	2
3	W 41-1688-6	Support post short	2
4	350-6203	Screwdriver Torx M4	1
5	240-7124	Screwdriver small No.1	1
6	PW 29.43.300	Pin for flexure holders	2
7	W 41-1700	Assembly jig with slot	1
8	W 41-1701	Assembly jig with hole	1
9	W 41-1677	Assembly pin for assembly jig	4
10	W 41-1691	Centring spigot for flexures	2
11	W 41-1774	Centring spigot for coil, EP/ES	1
12	W 41-1675	Centring spigot for coil, XT/XB	1
13	W 41-1674	Coupling anchor	1
14	PN 1100-172	Screw cyl. Torx M4x12	2
15	W 41-1699	Alignment jig	1
16	W 41-1578-8	Connecting cable	1
17	W 41-1681	Centring bush	1

4.1.4 Further tools needed


Instruments:

1. Soldering iron
2. Tweezers
3. Lens
4. Voltmeter

Standard tools:

1. Box spanner Nos.4, 5.5 and 7
2. Flat-bladed screwdriver Nos.1, 2 and 4
3. Phillips screwdriver No.2
4. Hexagon socket wrench WS 2.5mm

4.1.5 Preparing the tools


1. Screw the 2 long [2] and 2 short [3] support posts on the assembly plate.
2. Preheat the soldering iron to $310 \pm 5^{\circ}\text{C}$ ($590 \pm 10^{\circ}\text{F}$).

4.2 Opening a balance

4.2.1 General procedure for opening an XT/XB balance

1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 cylindrical screws on the housing's top using a No. 4 screwdriver.
4. Slightly raise the housing's top.
5. Release the membrane keypad ribbon cable from the main board.

4.2.2 General procedure for opening an IP65 balance


1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 screws IP56 with a No. 1 screwdriver and remove the clamping ring.
4. Unscrew the 3 cylindrical screws on the housing's top and the 2 cylindrical screws on the housing's bottom using a No. 4 screwdriver.
5. Slightly raise the housing's top (do not damage the sealing).
6. Release the membrane keypad ribbon cable from the main board.

4.3 Removing the reference weight device (SCS only)

1. Unplug the servo-motor (4/p.28) from the connection board (12/p.6), remove the 3 screws (8/p.28) and lift the complete device including its reference weight out of the housing's bottom.

4.4 Removing the weighing cell from the housing's bottom

1. Unplug the floating sensor plug (15/p.12) on the main board using a small screwdriver.
2. Turn over the housing's bottom and unscrew the 3 nuts (26/p.12) on the underside of the housing.
3. Hold the weighing cell only on the chassis and lift it out of the housing's bottom.
4. The chassis can now be placed on the assembly plate for easier handling.
5. To protect the weigh cell, fix the 2 assembly jigs [7] + [8] with the 4 assembly pins [9].

4.5 Dismantling the weighing cell

4.5.1 Removing the coupling

1. Unscrew the 2 screws (25/p.12), remove the coupling (3/p.12) incl. its bases (9/p.12) and spring washers (28/p.12) → Do not separate the screw components.
2. Screw the coupling anchor [13] tightly on with the 2 screws [14] in place of the coupling.

4.5.2 Removing the 2 flexure holders

1. Fix the upper flexure holder (2/p.12) into the chassis, using the 2 pins for flexure holders [6].
2. Unscrew the 4 screws (24/p.12) of the upper flexure holder and remove them together with its 4 spring washers (28/p.12) → 2 on the chassis, 2 on the support piece.
3. Remove the 2 pins with the upper flexure holder together with its 4 spring washers and flexures, then remove the pins from the holder.
4. Fix the lower flexure holder (2/p.12) into the chassis with the two pins for flexure holders [6].
5. Unscrew the 4 screws (24/p.12) from the lower flexure holder and remove them together with the 4 spring washers → 2 on the chassis, 2 on the support piece.
6. Remove the 2 pins with the lower flexure holder together with its 4 spring washers and flexures, then remove the pins from the holder.

4.5.3 Removing the coil

1. Unsolder the 2 connecting wires (1/p.12) carefully from the connecting board (17/p.12) using the tweezers.
2. Unscrew the 3 screws (23/p.12) from the magnet cover (8/p.12), remove the cover together with the transport safety device (4/p.12) and the spring washer (28/p.12).
3. Unscrew the 2 brass screws (21/p.12) from the coil (12/p.12) on the balance arm (19/p.12) using a No.2 screwdriver, remove the coil together with the washer (27/p.12).
4. Remove the coil carefully from the pot (34/p.12).

4.5.4 Removing the support piece

1. Unscrew the 3 screws (23/p.12) of the inner cone (12/p.12) and remove it.
2. Unscrew the upper screw (balance arm side) of the coupling anchor [13] and pull off the support piece together with the coupling anchor over the assembly pins [9].
3. Untighten the lower screw of the coupling anchor for roughly 2 turns. → The support piece and the coupling anchor are now loosely held together.

4.5.5 Removing the balance arm


1. Unscrew the brass screw (30/p.12) from the stopper (29/p.12), turn the stopper inwards ca. 90° and gently re-tighten the brass screw.
2. Unscrew the 4 screws (25/p.12) from the vertical flexures (6/p.12) and remove them together with its bases (11/p.12) and spring washers (28/p.12). **Note:** Since the lug of the floating sensor is firmly attached to the balance arm, the floating sensor must be unscrewed to dismantle the balance arm and gently tilted back:
3. Unscrew the 2 screws (20/p.12) from the floating sensor plate (5/p.12) and remove the plate.
4. Unscrew the screw (20/p.12) from the connection board (17/p.12).
5. Unscrew the 3 screws (22/p.12) of the floating sensor and gently tilt the floating sensor to the rear.
6. Remove the two inner assembly pins [9] → unscrew from the balance arm.
7. Remove the two outer assembly pins [9] → unscrew from the chassis.
8. Remove the two assembly jigs [7] + [8].
9. Carefully tilt the balance arm (19/p.12) over the pot and the floating sensor and take it out.
10. Further dismantling of the weighing cell is unnecessary.

4.6 Installing new horizontal flexures

4.6.1 Removing the old horizontal flexures

1. Remove all screws from the horizontal flexure holders and remove the old flexures together with their spring washers.


4.6.2 Installing new horizontal flexures


1. Place a flexure holder on the assembly plate [1] .
2. Select 4 new horizontal flexures (7/p.12) → Place them carefully on the assembly plate → their slot has to be always on the flexure holder.
3. Place the 4 spring washers (28/p.12) on the flexures and screw them tightly to the holder with 4 screws (22/p.12).
4. Take the complete flexure holder carefully off the assembly plate → Do not touch the flexures!
5. Prepare the other flexure holder the same way.

4.7 Cleaning


4.7.1 Cleaning the pot

1. Wrap a strips a double-sided adhesive tape around a flat stick.

2. Dab with the stick on the ground of the pot.

3. Prepare a thin and round stick the same way and spin it through the gap between the pot's wall and the magnet.
4. Dab over the top of the magnet pot with a piece of adhesive tape.

4.7.2 Cleaning the coil


1. Using a lint-free cloth, remove any contamination from the coil → Do not clean the coil with adhesive tape; risk of injury.

4.7.3 Checking the cleanliness again

1. Inspect the annular gap and the coil for residual contamination and clean if necessary.

4.8 Assembling the weighing cell

4.8.1 Assembling the balance arm

1. To facilitate the assembly, the chassis can be placed on the assembly plate → Pay attention to the floating sensor!
2. Place the assembly jigs [7] + [8] with 2 assembly pins [9] on the chassis, position them horizontally using the alignment jig [15] and tighten the two assembly pins [9] with a No. 4 box spanner.
3. Carefully introduce the balance arm (19/p.12) over the pot (34/p.12) and past it to the floating sensor.
4. Assemble the balance arm with 2 assembly pins [9] loosely behind the assembly jigs. The balance arm must remain moveable.
5. The right inner assembly pin of the assembly jig [8] can now be firmly tightened.
6. The left inner assembly pin of the assembly jig [7] can now be firmly tightened.


4.8.2 Assembling the coil, the magnet cover and the transport safety device

1. Place the coil onto the balance arm, centre it with the centring spigot [12] and fix it tightly with 2 brass screws (21/p.12) and 2 washers (27/p.12).
2. The centring spigot must be smoothly removable.
3. Pre-fix the magnet cover (8/p.12) with its left and right screw (23/p.12) .
4. Pre-assemble the transport safety device (4/p.12) with the third screw (23/p.12) and spring washer (28/p.12) → tighten gently.
5. The centring bush [17] shall be introduced between the balance arm-pilot and the transport safety device to assure constant annular gap. Tighten the screw (23/p.12) and remove the centring bush.

4.8.3 Assembling the floating sensor and the stopper


1. Press the floating sensor carefully into the cutout of the chassis in the direction of the pot and secure it to the chassis with the 3 screws (22/p.12).
2. Screw firmly the connection board (17/p.12) onto the chassis with the screw (20/p.12).
3. Fix the floating sensor plate (5/p.12) with the 2 screws (20/p.12) to the chassis.
4. Solder the two connection wires (1/p.12) parallelly and strain-freely to the connection board (17/p.12).
5. Turn back the stopper (29/p.12) with its brass screw (30/p.12) by 90° and tighten it.
6. The stopper has to be free of play and lie centrally in the symmetry screw (32/p.12) → if necessary, adjust the symmetry screw.

4.8.4 Assembling new vertical flexures


1. Select 2 new vertical flexures (6/p.12) see "Balance groups and spare parts – Waagengruppen und Ersatzteile" on page 4.
2. Bring the flexures into position → its slots have to be on the support piece side (below).
3. Centre the flexures in the hole on the chassis side (above) with a centring spigot [10].
4. Centre the flexures in the slot on the support piece side (below) with the second centring spigot [10].
5. Tighten the first (upper) centring spigot and then remove the second (lower) centring spigot. → Do not touch the flexures!
6. Fasten the flexures on the balance arm side with their bases (11/p.12), spring washers (28/p.12) and screws (25/p.12).
7. Remove the second (upper) centring spigot and fasten the vertical flexures on the chassis side with their bases (11/p.12), spring washers (28/p.12) and screws (25/p.12).

4.8.5 Assembling the support piece and the inner cone


1. Introduce the support piece with the loosely held coupling anchor over the two inner assembly pins [9].
2. It must be possible to push the support piece smoothly over the two inner assembly pins. In case the support piece can not be pushed smoothly running over the two assembly pins, unscrew the pin [9] and tighten it again.
3. Locate the support piece and the balance arm using the coupling anchor [13] and the two screws [14] → The support piece and balance arm are fixed now.
4. Assemble the inner cone (12/p.12) to the support piece with the 3 screws (23/p.12).

4.8.6 Assembling the flexure holders

1. Introduce the two pins for flexure holders [6] in the first pre-assembled holder (→ see "Installing new horizontal flexures" on page 59) the way that the shoulder of the pins is opposite the screw heads on the holder.
2. Bring the two pins with the holder into position in the chassis.
3. Lower the holder carefully over the two pins on the chassis, until the flexures just touch the chassis and the support piece.
4. Fix the holder with the 4 screws (24/p.12) and the 4 spring washers (28/p.12).
5. Push the two pins through, up to the installed horizontal flexures holder.
6. Introduce on the opposite side the second pre-assembled flexure holder (→ see "Installing new horizontal flexures" on page 59) with the screw heads inwards over the two pins and lower onto the chassis, until the flexures just touch the chassis and the support piece.
7. Secure the holder with the 4 screws (24/p.12) and the 4 spring washers (28/p.12).
8. It must be possible to pull out the two pins smoothly from the flexure holders.

4.8.7 Assembling a new coupling

1. Select a new coupling (3/p.12) see "Balance groups and spare parts – Waagengruppen und Ersatzteile" on page 4.
2. Remove the coupling anchor [13].
3. Bring the coupling into position → the slot is on the support piece side.
4. Centre the coupling first on the balance arm side, then on the support piece side, using a centring spigot [10].
5. Tighten the first (upper) centring spigot, then carefully remove the second (lower) centring spigot. Do not touch the coupling and fasten it on the support piece side with the base (11/p.12), the spring washer (28/p.12) and the screw (25/p.12) .
6. Remove the second (upper) centring spigot and fasten the coupling on the balance arm side with the base (11/p.12), the spring washer (28/p.12) and the screw (25/p.12).

4.8.8 Dismantling the assembly jigs

1. Remove the 2 inner assembly pins [9] carefully.
2. Remove the 2 outer assembly pins [9] carefully.
3. Take out the two assembly jigs [7] + [8] by sliding them sideways.
4. Remove the two assembly pins for flexure holders [6] carefully.

4.9 Installing the weighing cell into the housing's bottom

1. Carry out a further movement check with the lens on the balance arm → the balance arm must not touch the transport safety device (4/p.12)!
2. Place the weighing cell carefully into the housing's bottom → hold the cell only on the chassis.
3. Tilt the housing's bottom to the side and screw the weighing cell securely from below with the 3 nuts (26/p.12).
4. Connect the floating sensor plug (15/p.12) to the main board (11/p.6).

4.10 Installing the reference weight device into the housing's bottom (SCS only)

1. Introduce the device including its reference weight. **Caution:** Do not touch the weighing cell!
2. Align the device in the housing's bottom and secure with the 3 screws (8/p.28).
3. Connect the servo motor (4/p.28) to the connection board.

5. Adjustments

5.1 Before you get started


5.1.1 Table of specifications

1. In carrying out the following adjustments, it is necessary to know the specifications of the balance in service, which you may download from our website:
<http://www.precisa.com>
The following instructions refer often to that specifications, have them always at disposal!

5.1.2 Placing loads

1. Place all loads always carefully and, unless otherwise asked, in the center of the weighing pan!

5.2 Setting the balance to the factory mode (the balance is open)

1. Remove the jumper J1 from the main board.
2. Level the balance with the levelling screws and connect it to the mains with the power adapter → The balance starts up in the factory mode.

5.3 Adjusting the symmetry

1. Ensure the balance is in the factory mode → see "Setting the balance to the factory mode (the balance is open)" on page 68.
2. Connect the connecting cable [16] to plug S3 (see its position in the figure of the chapter 'Setting the balance to the factory mode (the balance is open)' on page 68) on the main board.
3. Connect the voltmeter to the connecting cable.
4. The voltmeter must show approximately the **same voltage** in both end positions of the balance arm (once positive, once negative), if it doesn't, adjust the position of the balance arm with the symmetry screw.
5. Detach the connecting cable and voltmeter from the balance.

5.4 Checking the pre-load


1. Ensure the balance is in the factory mode → see "Setting the balance to the factory mode (the balance is open)" on page 68.
2. The pre-load error ERROR 60 (below minimal converter range) must not appear with the weighing pan in position without load.
3. The pre-load error ERROR 61 (above maximal converter range) must not appear with the weighing pan in position with full load.
4. If such an error is displayed, the symmetry must be re-adjusted and the transport safety device adjusted from the start → see "Assembling the coil, the magnet cover and the transport safety device" on page 45.

5.5 Adjusting the corner load

5.5.1 Determining the corner load values



1. Ensure the balance is in factory mode → see "Setting the balance to the factory mode (the balance is open)" on page 68.
2. Place the upper part of the housing on the right alongside the lower part of the housing and connect the membrane keyboard to S4 (see its position in the figure of the chapter 'Setting the balance to the factory mode (the balance is open)' on page 68) on the main board.
3. Put the pan holder and the weighing pan on the balance.
4. Have, according to the specifications, a suitable corner load reference weight at disposal. For its placement positions refer to the figure of chapter 'Adjusting the corner load (example)' on page 71.
5. Place the weight on C, tare the balance, shift the weight to A and calculate the difference: $\Delta AC = A - C = A - 0 = A$ [d]; (C = tared = 0); **Example:** $\Delta AC = A = +12$ d.
6. Place the weight on D, tare the balance, shift the weight to B and calculate the difference: $\Delta BD = B - D = B - 0 = B$ [d]; (D = tared = 0); **Example:** $\Delta BD = B = -7$ d.

Note: Skip the following chapter 'Adjusting the corner load (example)' on page 71 if the determined corner load values stay within the tolerated range.

5.5.2 Adjusting the corner load (example)

1. Start adjusting the corner load on the diagonal with the larger difference.
2. The difference of diagonal A-C is positive (+12 d). Turning the corner load screw 1 clockwise reduces the value of A → the positive difference A-C becomes absolutely smaller.
3. The difference of diagonal B-D is negative (-7 d). Turning the corner load screw 2 anticlockwise raises the value of B → the negative difference B-D becomes absolutely smaller.
4. Continue determining the corner load error (as described in the previous chapter) and reducing it until both diagonal differences stay within the tolerated range.

5.6 Checking the lowering process of the internal reference weight (SCS only)

5.6.1 Checking the lowering process of the internal reference weight (XT only)

1. Ensure the balance is in factory mode → see "Setting the balance to the factory mode (the balance is open)" on page 68.
2. Place the upper part of the housing on the right alongside the lower part of the housing and connect the membrane keyboard to S4 (see its position in the figure of the chapter 'Setting the balance to the factory mode (the balance is open)' on page 68) on the main board.
3. Hold  until "TESTPROG 1" appears.
4. Lower the weight with **{DOWN}** and lift it with **{UP}**.
5. If the weight doesn't run cleanly into its holder, re-align the reference weight device in the housing's bottom.


5.6.2 Checking the lowering process of the internal reference weight (XB only)

1. Ensure the balance is in factory mode → see "Setting the balance to the factory mode (the balance is open)" on page 68.
2. Place the upper part of the housing on the right alongside the lower part of the housing and connect the membrane keyboard to S4 (see its position in the figure of the chapter 'Setting the balance to the factory mode (the balance is open)' on page 68) on the main board.
3. Hold **MODE** until "TESTPROG 1" appears.
4. Lower the weight with **{DOWN}** (hold **REF** until „FKT-KEY 4“ appears) and lift it with **{UP}** (hold **REF** until „FKT-KEY 3“ appears).
5. If the weight doesn't run cleanly into its holder, re-align the reference weight device in the housing's bottom.

5.7 Closing the balance

1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the weighing pan and the pan holder.
3. Ensure that the jumper J1 is set and the membrane keypad is connected by S4 to the main board (see their positions in the figure of the chapter 'Setting the balance to the factory mode (the balance is open)' on page 68).
4. Ensure that both, the display and the glass of the housing's top are clean.
5. Place the housing's top on its bottom and fasten together with the 3 screws (1x9/p.8 and 2x10/p.8).
6. **IP65 only:** Fix the housing with the two additional screws (10/p.11). Fix the membrane carefully with the clamping ring (4/p.11) and its 3 screws (2/p.11) without causing any wrinkles on the membrane.
7. Put the weighing pan and the pan holder and, if present, the windshield on the balance.

5.8 Setting the balance to the service mode (the balance is closed)

5.8.1 Setting the balance to the service mode (XT only)

1. Level the balance with the levelling screws and connect it to the mains with the power adapter.
2. Press the service switch on the rear of the housing's bottom with the small screwdriver [5] → The display flashes once.
3. Press  → „CODE 0“ appears.
4. Press  → „CODE 0000“ appears.
5. Enter the service access code 5236 with ,  and  and confirm with  → The balance is now in the service mode.

5.8.2 Setting the balance to the service mode (XB only)

1. Level the balance with the levelling screws and connect it to the mains with the power adapter.
2. Press the service switch on the rear of the housing's bottom with the small screwdriver [5] → The display flashes once.
3. Press **MODE** twice → „CODE 0000“ appears.
4. Pressing **MODE** changes the value of the figure, holding changes to the next figure.
5. Enter the service access code 5236 and confirm with **PRINT** → The balance is now in the service mode.

5.9 Adjusting the linearity

IMPORTANT: Adjusting the linearity erases all S-correction factors, too. → This chapter has therefore to be processed entirely!

5.9.1 Adjusting the linearity (XT only)

1. Ensure the balance has reached its operating temperature and is in service mode → see "Setting the balance to the service mode (XT only)" on page 73.
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Hold  until „LINEARITY“ appears.
4. Have the weighing pan empty and press **{0}** → the balance catches the zero point and the display flashes.
5. Wait until the flashing stops, place the half load and press **{1/2}** → the balance catches the half load value and the display flashes.
6. Wait until the flashing stops, place the full load and press **{1/1}** → the balance catches the full load value and the display flashes.
7. Wait until the flashing stops and press **{SET}** → the balance sets the linearisation factors and displays the full weight value.
8. Empty the weighing pan, tare and check if the linearity stays within the tolerated range. If it doesn't, process chapter 'Adjusting the S-correction' on page 77.
9. **With SCS:** Process chapter 'Adjusting the internal reference weight (SCS only)' on page 79.
Without SCS: Switch the balance off and on (by pressing **ON/OFF** twice) → the balance quits the service mode.

5.9.2 Adjusting the linearity (XB only)

1. Ensure the balance has reached its operating temperature and is in service mode → see "Setting the balance to the service mode (XT only)" on page 73.
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Hold **MODE** until „LINEARITY“ appears.
4. Have the weighing pan empty and press **{0}** (hold **REF** until „FKT-KEY 1“ appears) → the balance catches the zero point and the display flashes.
5. Wait until the flashing stops, place the half load and press **{1/2}** (hold **REF** until „FKT-KEY 2“ appears) → the balance catches the half load value and the display flashes.
6. Wait until the flashing stops, place the full load and press **{1/1}** (hold **REF** until „FKT-KEY 3“ appears) → the balance catches the full load value and the display flashes.
7. Wait until the flashing stops and press **{SET}** (hold **REF** until „FKT-KEY 4“ appears) → the balance sets the linearisation factors and displays the full weight value.
8. Empty the weighing pan, tare and check if the linearity stays within the tolerated range. If it doesn't, process chapter 'Adjusting the S-correction' on page 77.
9. **With SCS:** Process chapter 'Adjusting the internal reference weight (SCS only)' on page 79.
Without SCS: Switch the balance off and on (by pressing **ON/OFF** twice) → the balance quits the service mode.

5.9.3 Checking the linearity

1. Ensure the balance has reached its operating temperature.
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Tare the balance and check if the linearity stays at half, full, quarter and three quarter load within the tolerated range. If it doesn't at half or full load, adjust the linearity again. If it doesn't at quarter or three quarter load, process chapter 'Adjusting the S-correction' on page 77.

5.10 Checking the adjustments of the closed balance

5.10.1 Checking the repeatability (with SCS)

1. Ensure the balance has reached its operating temperature.
2. Have, according to the specifications, a suitable full load reference weight at disposal.
3. Ensure by the configuration menu (see operating instructions) the calibration/adjustment mode is set to „INTERNAL“.
4. Adjust the balance automatically by holding **TARE** until „CALIBRATION“/„ADJUSTMENT“ appears.
5. Tare the balance and check if the repeatability stays within the tolerated range. If it doesn't, process chapter 'Adjusting the internal reference weight (SCS only)' on page 79 and process this chapter completely again.

5.10.2 Checking the repeatability (without SCS)

1. Ensure the balance has reached its operating temperature.
2. Have, according to the specifications, a suitable full load reference weight at disposal.
3. Ensure by the configuration menu (see operating instructions) the calibration/adjustment mode is set to „EXTERNAL“.
4. Adjust the balance manually with the ready full load weight by holding **TARE** until „CALIBRATION“/„ADJUSTMENT“ appears and following the displayed instructions, see also the operating instructions.
5. Tare the balance and check if the repeatability stays within the tolerated range. If it doesn't, re-adjust the balance and process this chapter completely again.

5.10.3 Checking the corner load


1. Have, according to the specifications, a suitable corner load reference weight at disposal. For its placement positions refer to the figure of chapter 'Adjusting the corner load (example)' on page 71.
2. Place the weight on C, tare the balance, shift the weight to A and calculate the difference: $\Delta AC = A - C = A - 0 = A$ [d]; (C = tared = 0).
3. Place the weight on D, tare the balance, shift the weight to B and calculate the difference: $\Delta BD = B - D = B - 0 = B$ [d]; (D = tared = 0).
4. Check if the determined corner load values stay within the tolerated range. If it doesn't, re-adjust it as described in chapter 'Adjusting the corner load (example)' on page 71 and process this chapter completely again.

5.10.4 Checking the linearity

1. Process chapter 'Checking the linearity' on page 75.

5.11 Adjusting the S-correction

5.11.1 Measuring the S-curve

1. Ensure the balance has reached its operating temperature and is in weighing mode (turn it off and on by pressing **ON/OFF** twice).
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 1x half load, 1x quarter load.
3. Tare the balance, place the quarter load and note the displayed value.
4. Place the three quarter load and note the displayed value.
5. Calculate the differences:
 $\Delta a = (\text{displayed quarter load}) - (\text{placed quarter load})$
 $\Delta b = (\text{displayed three quarter load}) - (\text{placed three quarter load})$

Example:

$$\Delta a = 50.0004g - 50.0000g = 0.0004g; \Delta b = 149.9998g - 150.0000g = -0.0002g$$

5.11.2 Calculating the S-correction factors

1. Indicate the differences in digits and change their signs:
 $Sa = -\Delta a [d]$
 $Sb = -\Delta b [d]$

Example:

$$\Delta a = 0.0004g = +4 d \rightarrow Sa = -\Delta a = -4 d$$

$$\Delta b = -0.0002g = -2 d \rightarrow Sb = -\Delta b = +2 d$$

5.11.3 Entering the S-correction factors (XT only)

1. Ensure the balance is in service mode → see "Setting the balance to the service mode (XT only)" on page 73.
2. Hold  until "SCORRECTION", followed by „SET SECTOR-A +0" appears.
3. Enter the S-correction factor Sa with  and  and set it with **{SET}** (press ).
4. Switch to sector b with **{SECTOR}** (press  → „SET SECTOR-B +0" appears.
5. Enter the S-correction factor Sb with  and  and set it with **{SET}** (press ).
6. Hold  until „BALANCING" appears.
7. Process chapter 'Checking the adjustments of the closed balance" on page 76.

5.11.4 Entering the S-correction factors (XB only)

1. Ensure the balance is in service mode → see "Setting the balance to the service mode (XB only)" on page 73.
2. Hold **MODE** until "SCORRECTION", followed by „SET SECTOR-A +0" appears.
3. Enter the S-correction factor Sa
with **{UP}** (hold **REF** until „FKT-KEY 3" appears)
and **{DOWN}** (hold **REF** until „FKT-KEY 4" appears)
and set it with **{SET}** (hold **REF** until „FKT-KEY 1" appears).
4. Switch to sector b with **{SECTOR}** (hold **REF** until „FKT-KEY 2" appears) → „SET SECTOR-B +0" appears.
5. Enter the S-correction factor Sb
with **{UP}** (hold **REF** until „FKT-KEY 3" appears)
and **{DOWN}** (hold **REF** until „FKT-KEY 4" appears)
and set it with **{SET}** (hold **REF** until „FKT-KEY 1" appears).
6. Hold **MODE** until "BALANCING" appears.
7. Process chapter 'Checking the adjustments of the closed balance" on page 76.

5.12 Adjusting the internal reference weight (SCS only)

It is possible to vary the full load reference weight within particular step sizes which may be, according to your actual balance, 10, 50, 100 or 500g. This allows an adjustment whether with various reference weights at disposal or close to the usual operating range of the balance.

IMPORTANT: Ensure the reference weight lies within the range of $\pm 0.3d$ of its nominal value!

Example: A 4200C of a customer, who usually weighs loads of 2.8kg, asks for 2000g. Its reference weight step size is 500g. A reference weight of 3000g (2000g + 2 steps of 500g) is the closest to the usual operating range. → Place such a reference weight if its mass lies within the range of 2999.997g .. 3000.003g (maximum deviation = $\pm 0.3d = \pm 0.003g$ at $d=0.01g$, see specifications).

1. Ensure the balance has reached its operating temperature and is in service mode:
XT: see "Setting the balance to the service mode (XT only)" on page 73.
XB: see "Setting the balance to the service mode (XB only)" on page 73.
2. Have, according to the specifications, a suitable full load reference weight at disposal.
3. Hold **TARE** until „CALIBRATION“/„ADJUSTMENT“ appears → the balance starts carrying out an internal adjustment without setting the obtained correction factors, starts thereupon with the external adjustment and asks, as soon the zero point is measured, for the reference weight. Its value is displayed flashing.
4. Place the asked weight, or better, the one placed at disposal → the balance catches its value and the display flashes. As soon as the flashing stops, the balance asks for the zero point. This is displayed flashing: „-- 0000 g“.
5. Remove the weight → the balance continues with the external adjustment, determines newly the correction factors for the internal reference weight, carries out immediately an internal adjustment and finally turns itself off and on, in order to quit the service mode.

5.13 Error messages

5.13.1 Operating error; Error disappears with correct use or application

Code	Cause	Remedy
1	Invalid Firmware	Load right firmware
3	Load too small	Use higher load
4	Wrong anti theft code	Enter right code
5	7 times wrong anti theft code	Enter unlock code → contact service
8	Zero point out of range	
9.0	Bad adjustment factor	Repeat adjustment
9.1	Zero point deviation exceeded	Repeat adjustment
9.2	Deviation of int. ref. wgt. exceeded	Repeat adjustment
9.3	Bad adjustment factor	Repeat adjustment
30	No internal reference weight	Acknowledge with TARE Adjust internal reference weight
31	No external reference weight	Acknowledge with TARE Define external reference weight
47	Too much values	Record less values
48	Not the pristine unit	Set to pristine unit
49	Too less values	Record more values
53	Values out of the +/-50% range	Place load within the range
60	Converter range fallen below	
61	Converter range exceeded	

5.13.2 Fatal error; Balance stays still, error remediable in service mode

Code	Cause	Remedy
10	No linearisation values	Repeat linearisation
11	No adjustment factor	Repeat adjustment
21	No SCS compens. values, Tmp.	Contact service
22	No compens. values, Tmp.	Contact service
23	No SCS compens. values, Tmp., Nonlin.	Contact service

5.13.3 Hardware error; Balance stays still

Code	Cause	Remedy
14	Compatibility conflict	Contact service
16	Broken internal RAM	Exchange controller
17	Broken external RAM	Exchange RAM
18.1	Wrong firmware checksum	Load right firmware Exchange flash memory
18.2	No approved firmware in approved balance	Load right firmware
18.8	Empty flash memory	Load right firmware Exchange flash memory
18.9	Wrong checksum of controller firmware	Exchange controller
19	Compatibility conflict	Contact service
20	Erroneous temperature measurement	Check whether floating sensor connected

5.14 Discussion with the manufacturer

IMPORTANT: Have at disposal the following information in discussions with the manufacturer:

1. Balance model: *e.g. XT 220A*
2. Order number: *e.g. 320-9201 (Label on the rear of the balance)*
3. Serial number: *e.g. 2405-037 or 2805037 (Label on the rear of the balance)*
4. Sales number: *e.g. N 43210 (Label on the rear of the balance)*
5. Software ID *e.g. N10-0400.P09 (appears on the display when switching on the balance)*
6. A precise description of the failure
7. Weighing results

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